



PRODUCT DATA SHEET

Miniplast Ein-Shemer ACS Ltd. Serological Pipette, 1 ml

Miniplast Cat. No. 65001-01, 65001-02

1. Description/Specification

- 1.1 Description: One-piece construction, crystal clear PS serological pipette, graduated, edged color-coded band, filter, 1 ml, packed in paper/plastic wrap for easy opening
- 1.2 Dimensions: Length: 270 mm (± 2 mm)
- 1.3 Volume: 1 ml, accuracy $\pm 2\%$ at full volume
 - 1.3.1 Negative graduations: up to -0.17 ml
- 1.4 Material
 - 1.4.1 Pipette: crystal clear PS, free from heavy metals (EU directive 2000/53/EC and its amendments)
 - 1.4.2 Filter: Bonded polyolefin fibers
- 1.5 Level of transparency of pipette: crystal clear
- 1.6 Color
 - 1.6.1 Filter: white
 - 1.6.2 Identifying band: yellow
- 1.7 Sterilization: By Gamma ray
- 1.8 Quality control: Testing of various characteristics in accordance with the valid specification
- 1.9 Other information: For single use only; not to be pipetted by mouth; inhalation/ingestion to be avoided

2. Features

- 2.1 Non-cytotoxic (ISO 10993-5)
- 2.2 Non-pyrogenic (US Pharmacopeia 32, National formulary 27 [USP], General chapter 85)
- 2.3 Non-hemolytic (ISO 10993-4)
- 2.4 DNase-, RNase-free
- 2.5 Autoclavability: No
- 2.6 Storing conditions: Room temperature
- 2.7 Chemical resistance: See website:
www.miniplast.co.il/about.php/downloads/quality_assurance
- 2.8 Shelf life: 5 years after month of production

3.

Packaging	65001-01	65001-02
3.1 Units/Package	1	25
3.2 Units/Case	500	1000
3.3 Units/Master Carton	1500	3000