



PRODUCT DATA SHEET

Miniplast Ein-Shemer ACS Ltd.
Disposable Cup, 120 ml
Miniplast Cat. No. 40120-xx
Valid for Cat. Nos.: as per items in chart below

1. Description/Specification

- 1.1 Description: Disposable PP cup, 120 ml
- 1.2 Dimensions:
 - 1.2.1 Cup: Ø 56.8 mm at rim, Ø 48.6 at base; H 72.5 mm
 - 1.2.2 Cover: Ø 61.0 mm
- 1.3 Volume: 120 ml
- 1.4 Material:
 - 1.4.1 Cup: Polypropylene (PP)
 - 1.4.1.1 Complies with requirement of Directive 2002/72/EC and requirement of FDA Regulation 21 CFR 177.1640
 - 1.4.1.2 Free from heavy metals according to ELV directive 2000/53/EC and its following amendments
 - 1.4.2 Cover: High Density Polyethylene (HDPE)
- 1.5 Color:
 - 1.5.1 Cup: translucent
 - 1.5.2 Cover: red
 - 1.5.3.1 Other colors available with minimum ordering quantities
- 1.6 Degree of cleanliness: aseptic (machine sterile) when packaged in cellophane bag (1 or 5 units per bag)
- 1.7 Quality control: Testing of various characteristics in accordance with the valid specification
- 1.8 Other information: For single use only

2. Features

- 2.1 Basic features: Hydrophobic, embossed graduation and writing area, screw cover, leak-resistant design
- 2.2 Storing conditions: Room temperature
- 2.3 Chemical resistance: See website:
www.miniplast.co.il/about.php/downloads/quality_assurance

3. Packaging	40120-01	40120-40	40120-41	40120-43	40120-44	4112041
Units/bag	-	1	5	5	-	-
Units/carton	1200	250	300	480	700+700	1200 covers
Cover	no	yes	yes	yes	yes	only
Cover assembled	no	yes	yes	yes	no	no