



PRODUCT DATA SHEET

Miniplast Ein-Shemer ACS Ltd.
Disposable Cup and Cover
Miniplast Cat. No. 40200-xx
Valid for Cat. Nos.: 40200-01, 40200-02, 40200-03
and 41000-xx (Cover)
Valid for 41000-02, 41000-03, 41000-20

1. Description/Specification

- 1.1 Description: Disposable PP cup, 200 ml
- 1.2 Dimensions:
 - 1.2.1 Cup: Ø 69.6 mm; H 87.8 mm
 - 1.2.2 Cover: Ø 70.0 mm
- 1.3 Volume: 200 ml
- 1.4 Material:
 - 1.4.1 Cup: Polypropylene (PP)
 - 1.4.1.1 Complies with requirement of Directive 2002/72/EC and requirement of FDA Regulation 21 CFR 177.1640
 - 1.4.1.2 Free from heavy metals according to ELV directive 2000/53/EC and its following amendments
 - 1.4.2 Cover: High Density Polyethylene (HDPE)
 - 1.4.2.1 Complies with requirements of Commission Regulation (EC) no. 975/2009, amending directive 2012/72 EC; for translucent cover only: free from heavy metals as per cup
- 1.5 Color:
 - 1.5.1 Cup: translucent
 - 1.5.2 Cover: translucent or green (as per cat. no. details)
- 1.6 Degree of cleanliness: aseptic (40200-02 only)
- 1.7 Quality control: Testing of various characteristics in accordance with the valid specification
- 1.8 Other information: For single use only

2. Features

- 2.1 Basic features: Hydrophobic, graduation and writing area, push-down cover, leak-resistant design
- 2.2 Storing conditions: Room temperature
- 2.3 Chemical resistance: See website:
www.miniplast.co.il/about.php/downloads/quality_assurance

3. Packaging	40200-01	40200-02 (aseptic)	40200-03
Units per inner bag	1000	8	1000
Units per carton	1000	280	1000
With attached, green cover	no	yes	no
With separate, clear cover	no	no	yes

3a. Packaging of cover	41000-02	41000-03	41000-20
Units per inner bag	2500	100	2000
Units per carton	2500	2000	2000