

PRODUCT DATA SHEET

Miniplast Ein-Shemer ACS Ltd. Disposable PS Centrifuge Tube Miniplast Cat. No. 30050-xx

1. Description/Specification

- 1.1 Description: Disposable PS centrifuge tube, non-sterile
 - 1.1.1 Tube: Ø 29.5 mm; H 113.5 mm
 - 1.1.2 Tube with cap: Ø 34.7 mm; H 116.5 mm
- 1.2 Volume: 50 ml
- 1.3 Material:
 - 1.3.1 Tube: Polystyrene (PS), free from heavy metal
 - 1.3.2 Cap: High Density Polyethylene (HDPE)
- 1.4 Color:
 - 1.4.1 Tube: transparent (crystal clear)
 - 1.4.2 Printing on tube: white (graduation, writing area)
 - 1.4.3 Cap: red
 - 1.5.3.1 Other colors available with minimum ordering quantities
- 1.5 Sterilization: none
- 1.6 Quality control: Testing of various characteristics in accordance with the valid specification
- 1.7 Other information: Recommended for single use only

2. Features

- 2.1 Basic features: Hydrophobic, conical/skirted bottom, printed graduation and writing area, screw cap, leak-resistant design
- 2.2 Temperature range:
 - 2.2.2 Tube: +4° C to +80°C
 - 2.2.3 Cap: -80°C to +121°C
- 2.3 Centrifugation, max: 3080 RCF
- 2.4 Chemical resistance: See website:
www.miniplast.co.il/about.php/downloads/quality_assurance

3. Packaging: According to clients' requirements and minimum quantities