



## PRODUCT DATA SHEET

### Miniplast Ein-Shemer ACS Ltd. Disposable Petri Dish 90 mm Miniplast Cat. No. 20090-01, 25090-01, 20090-02

#### 1. Description/Specification

- 1.1 Description: Disposable PS Petri Dish with 3-vent lid
- 1.2 Dimensions: Ø 90 x 17 mm
- 1.3 Material: Polystyrene (PS) (complies with requirement of Directive 2002/72/EC and requirement of FDA Regulation 21 CFR 177.1640)
- 1.4 Color: Crystal clear
- 1.5 Degree of cleanliness: Aseptic (machine sterile)
- 1.6 Quality control: Testing of various characteristics in accordance with the valid specification
- 1.7 Other information:
  - 1.7.1 For single use only
  - 1.7.2 One-sector dish can be irradiated by gamma ray; Cat. No. 25090-01
  - 1.7.3 Available in two sectors; Cat. No. 29090-02

#### 2. Features

- 2.1 Basic features: Hydrophobic
- 2.2 Storing conditions: Room temperature
- 2.3 Chemical resistance: See website:  
[www.miniplast.co.il/about.php/downloads/quality\\_assurance](http://www.miniplast.co.il/about.php/downloads/quality_assurance)
- 2.4 Shelf life: 3 years after month of production (sterile), 1 year (aseptic)

#### 3.

Packaging	20090-01 Aseptic One sector	25090-01 Sterile One sector	20090-02 Aseptic Two Sectors
3.1 Units/Sleeve	20	20	20
3.2 Units/Carton	500	500	500