



PRODUCT DATA SHEET

Miniplast Ein-Shemer ACS Ltd.
Disposable Spreader
Miniplast Cat. No. 15004-xx
Valid for Cat. Nos.: 15004-01, 15004-05

1. Description/Specification

- 1.1 Description: Disposable PS spreader
- 1.2 Dimensions: 184 mm x 41.6 mm
- 1.4 Material: Polystyrene (PS)
- 1.5 Color: Green
- 1.6 Sterilization: By Gamma ray; SAL 10^{-3}
- 1.7 Quality control: Testing of various characteristics in accordance with the valid specification
- 1.8 Other information: For single use only

2. Features

- 2.1 Basic features: Hydrophobic; T-shaped
- 2.2 Storing conditions: Room temperature
- 2.3 Chemical resistance: See website:
www.miniplast.co.il/about.php/downloads/quality_assurance
- 2.4 Shelf life: 3 years after month of production

3.

Packaging	15004-01	15004-05
3.1 Units/Packet	1	5
3.2 Units/Case	80	200
3.3 Cases/Master Carton	10	10