

## PRODUCT DATA SHEET

**Miniplast Ein-Shemer ACS Ltd.**  
**QuadLoop**  
**Disposable Inoculation Loop**  
**Miniplast Cat. No. 15012-xx**

**Valid for Cat. Nos.: 15012-10, 15012-20, 15012-50**

### 1. Description/Specification

- 1.1 Description: Disposable PS inoculation loop 10 µl/needle
- 1.2 Dimensions: Length: 200 mm
- 1.3 Volume (wetting): 10 µl
- 1.4 Material: Polystyrene (PS)
- 1.5 Color: Orange
- 1.6 Sterilization: By Gamma ray; SAL 10<sup>-3</sup>
- 1.7 Quality control: Testing of various characteristics in accordance with the valid specification
- 1.8 Other information: For single use only

### 2. Features

- 2.1 Basic features: Hydrophobic; unique, four-sided handle arrangement
- 2.2 Storing conditions: Room temperature
- 2.3 Chemical resistance: See website:  
[www.miniplast.co.il/about.php/downloads/quality\\_assurance](http://www.miniplast.co.il/about.php/downloads/quality_assurance)
- 2.4 Shelf life: 3 years after month of production

### 3.

<b>Packaging</b>	<b>15012-10</b>	<b>15012-20</b>	<b>15012-50</b>
3.1 Units/Packet	10	20	50
3.2 Units/Case	720	1000	1000
3.3 Cases/Master Carton	10	10	10