



PRODUCT DATA SHEET

Miniplast Ein-Shemer ACS Ltd.
QuadLoop
Disposable Inoculation Loop
Miniplast Cat. No. 15002-xx

Valid for Cat. Nos.: 15002-10, 15002-20, 15002-50

1. Description/Specification

- 1.1 Description: Disposable PS inoculation loop 1 µl/10 µl
- 1.2 Dimensions: Length: 200 mm
- 1.3 Volume (wetting): 1 µl/10 µl
- 1.4 Material: Polystyrene (PS)
- 1.5 Color: Blue
- 1.6 Sterilization: By Gamma ray; SAL 10⁻³
- 1.7 Quality control: Testing of various characteristics in accordance with the valid specification
- 1.8 Other information: For single use only

2. Features

- 2.1 Basic features: Hydrophobic; unique, four-sided handle arrangement
- 2.2 Storing conditions: Room temperature
- 2.3 Chemical resistance: See website:
www.miniplast.co.il/about.php/downloads/quality_assurance
- 2.4 Shelf life: 3 years after month of production

3.

Packaging	15002-10	15002-20	15002-50
3.1 Units/Packet	10	20	50
3.2 Units/Case	720	1000	1000
3.3 Cases/Master Carton	10	10	10